




## ***Wireless Remote Programming Guide***

- 1. If you are only reprogramming the wireless remote skip this step! Remove the Philips head screw from the boom holding the remote lanyard to the boom. Once removed, re-install the Philips head screw with the new remote/lanyard to the boom.***



- 2. Unplug the negative battery cable (black wire) from the battery terminal.***



- 3. Make sure the control box key is in the ON position.**
- 4. While pressing and holding the up arrow on your wireless remote, plug in the negative battery cable (black wire) to the negative post on the battery.**



- 5. Within 5 seconds after re-installing the negative battery cable to the terminal, the lift should start to raise the chair/sling.**

**6. Make sure to lower/raise the chair, and rotate a full 360° to test the new wireless remote.**